



## **FB Plasmid Maxiprep Kit**

www.fairbiotech.com

For research use only

Sample: up to 200 ml bacterial cells

Yield: up to 850µg of plasmid

#### Introduction

The FB Plasmid Maxiprep Kit provides a fast, simple, and cost-effective plasmid maxiprep method for isolation of plasmid DNA from cultured bacterial cells. The FB Plasmid Maxiprep Kit is based on alkaline lysis of bacterial cells followed by binding of DNA onto the glass fiber matrix of the spin column in the presence of a high amount of salt.. Phenol extraction and ethanol precipitation are not required, and high-quality plasmid DNA is eluted with a small volume of Tris buffer (included in each kit) or water (pH is between 7.0 and 8.5). Plasmid DNA purified with FB Plasmid Maxiprep Kit is suitable for a variety of routine applications including restriction enzyme digestion, Sequencing, library screening, in vitro translation, transfection of robust cells, ligation and transformation. The entire procedure can be completed within 40-50 minutes.

#### **Kit Contents**

Catalog No.	XP02	XP10	XP20
Buffer M1	20 ml	85 ml	85 ml X 2
Buffer M2	20 ml	85 ml	125 ml , 40 ml
Buffer M3	25 ml	125 ml	125 ml X 2
Buffer W1	25 ml	105 ml	105 ml X 2
Buffer W2	6 ml	25 ml	25 ml X 2
(Add Ethanol)	(24 ml)	(100 ml)	(100ml) X 2
Buffer EL	5 ml	30 ml	50 ml
RNase A (50mg/ml)	Added	200µl	200µl X 2
XP Columns	2 pcs	10 pcs	20 pcs

## **Quality Control**

In accordance with FairBiotech's ISO-certified Total Quality Management System, the quality of the FB Plasmid Maxiprep Kit is tested on a lot-to-lot basis to ensure consistent product quality.

### **Additional requirements**

- \*Ethanol (96~100%)
- \*50 ml centrifuge tubes

#### **NOTE**

- ★ Add the provided RNase A solution to Buffer M1, mix, and store at 2-8°C.
- ★ Add ethanol (96–100%) to Buffer W2 before use (see bottle label for volume).
- ★ Check Buffers before use for salt precipitation. Redissolve any precipitate by warming to 37°C
- ★ Buffers M2, M3, and W1 contain irritants. Wear gloves when handling these buffers.

# **FB Plasmid Maxiprep Kit Protocol**

## **Step 1 Bacterial Cells Harvesting**

- Transfer 200 ml bacterial culture to a centrifuge tube.
- Centrifuge at 6,000 x g for 5 minute and discard the supernatant.

#### Step 2 Resuspend

Resuspend pelleted bacterial cells in 8 ml Buffer M1 (RNase A added)

## Step 3 Lysis

Add 8 ml Buffer M2 and mix thoroughly by inverting the tube 10 times (Do not vortex) and then stand at room temperature for 2 minutes or until the lysate is homologous.

## **Step 4 Neutralization**

Add 12 ml Buffer M3 and mix immediately and thoroughly by inverting the tube 10 times (Do not vortex). Centrifuge at 6,000x g for 15 minutes.

#### Step 5 Binding

♦ Place a XP Column in a 50 ml centrifuge tube. Apply 15 ml of the supernatant (from step 4) to the XP column by decanting

- or pipetting.
- Centrifuge at 6,000 x g for 3 minutes. Discard the flow-through and place the XP column back into the same 50 ml centrifuge tube. Transfer the remaining supernatant to the same XP Column.
- Centrifuge at 6,000 x g for 3 minutes. Discard the flow-through and place the XP column back into the same 50 ml centrifuge tube.

## Step 6 Wash

- ◆ Add 10 ml of **Buffer W1** into the **XP Column**. Centrifuge at 6,000 x g for 3 minutes. Discard the flow-through and place the **XP column** back into the same 50 ml centrifuge tube.
- Add 12 ml of Buffer W2 (Ethanol added) into the XP Column. Centrifuge at 6,000 x g for 3 minutes. Discard the flow-through and place the XP column back into the same 50 ml centrifuge tube.
- ◆ Centrifuge at 6,000 x g again for 3 minutes to remove residual **Buffer W2**.

#### Step 7 Elution

- ◆ To elute DNA, place the XP column in a new 50 ml centrifuge tube.
- Add 2 ml Buffer EL or water (pH is between 7.0 and 8.5) to the center of each XP column, let stand for 2 minutes, and centrifuge at 6,000 x g for 3 minutes.



## **Troubleshooting**

Problem	Cause	Solution	
Presence of RNA	RNA contamination	Prior to using Buffer M1, ensure RNase A is added.	
Plasmid bands was smeared on agarose gel	plasmid DNA degradation	Keep plasmid preparations on ice or frozen in order to avoid plasmid DNA degradation	
Presence of genomic DNA	Genomic DNA contamination	Do not overgrow bacterial cultures.  Do not incubate more than 5 min after adding Buffer M1.	
Low yields of DNA	Low plasmid copy number	Increase the culture volume. Change the culture medium.	
	96~100% ethanol not used	Add ethanol (96~100%) to the Buffer W2 before use.	
	Nuclease contamination	Check buffers for nuclease contamination and replace if necessary.  Use new glass- and plastic-wares, and wear gloves.	
	Column overloaded	Decrease the loading volume or lower the culture density.	
	SDS in the Buffer S2 precipitated	The SDS in Buffer M2 may precipitate with storage. If this happens, incubate the Buffer M2 at 30~40°C for 5 min and mix well.	
	Incorrect elution conditions	Ensure that Buffer E is added into the center of the MP Column.	
	Plasmid lost in the host E. coli	Prepare and use fresh culture.	



Inhibition of downstream	TE buffer used for DNA	Use the ethanol to precipitate the DNA, or repurify the DNA
enzymatic reactions	elution.	fragments and elute with the nuclease-free water.
	Presence of residual ethanol	Following the Wash Step, dry the MP Column with an
	in plasmid.	additional centrifugation step at 6,000 x g for 5 minutes.
DNA passed through in the	Column overloaded	Check the culture volume. If overgrown, add additional
flow-through or wash fraction		reaction buffer.
		Check the loading volume.
	Inappropriate salt or pH	Ensure that any buffer prepared in the laboratory was
	conditions in buffers	prepared according to instructions.
Plasmid DNA floats out of	Incomplete removal of the	Make sure that no residual ethanol remains in the
wells while running in	ethanol	membrane before eluting the plasmid DNA.
agarose gel		Re-centrifuge or vacuum again if necessary.

