



# MorreTaq High-Fidelity DNA Polymerase

## The High-Fidelity Polymerase For Longer DNA

MorreTaq High-Fidelity DNA polymerase are thermostable enzymes formulation specifically developed for synthesizes length up to 30kb and low error rate PCR product.

High-Fidelity polymerase synthesizes higher yields of product from genomic DNA, cDNA, bacterial cultures. It is contain a 2.5 hours half life at 96oC and easy to amplify PCR product at G-C rich and secondary structure.

### 10X reaction buffer

Containing 25mM MgCl<sub>2</sub>.

### Unit description

One unit of MorreTaq High-Fidelity DNA polymerase incorporates 10 nmol of dNTP into acid-insoluble material in 30 min at 74°C.

### Storage conditions

Storage at -20°C

### Template

MorreTaq High-Fidelity DNA polymerase is suitable for amplifying targets up to 15 kb from the following templates:

Genomic DNA: 10–200 ng

Plasmid DNA : 1–5 ng

cDNA : ~100 ng starting total RNA

Amplification of longer targets (up to 15 kb) may be possible, but may require more template and longer elongation times.

### Primers

Use 0.3 μM per primer as a general starting point.

For larger amounts of template (e.g., 200 ng genomic DNA), increasing the concentration up to 0.5 μM per primer may improve yield.

### Annealing Temperature

The annealing temperature is slightly higher than with typical PCR. The optimal annealing temperature should be ~2°C lower than the T<sub>m</sub> of the primers used. A range of 58–68°C is recommended.

### Extension Time

As little as 30 seconds per kb is suitable for most targets. Use up to 60 seconds per kb for maximum yield.

### PCR reaction mix

Component	Volume
MorreTaq High-Fidelity DNA polymerase	0.5-1ul
10X buffer	10 ul
10mM dNTP	2 ul
Primer1 (20 pmol)	2-4 ul
Primer2 (20 pmol)	2-4 ul
template	1-10 ul
ddH <sub>2</sub> O	Up to 100 ul
Total	100 ul

### PCR cycles

Step	Temperature	Time	Cycle
Initial denaturation	94-96°C	0.5-2 mins	1
Denaturation	94-96°C	0.2-2 mins	
Annealing	50-68°C	0.2-2 mins	15-30
Extension	68-75°C	1min/1kb	
Final extension	68-75°C	1-10 mins	1

### Step

After cycling, maintain the reaction at 4°C. Samples can be stored at -20°C until use.

### Ordering information

Cat.	Pack	Con.
MHF500	500U	5U/ul
MHF2500	2500U	5U/ul

For Research Use Only.

Not for use in diagnostic procedures.

About **MORREBIO**

Find more on <https://www.light-biotech.com/>